




Starkey | Technical Data

RIC

Receiver-In-Canal



Fitting Range



Genesis AI 24 | 20 | 16 | 12

Color Guide

Standard Colors



Patient Features

- IP Rating: 68
- Tinnitus Technology
- Wireless Connectivity
- CROS System
- Dual Radio (2.4 GHz + NFMI)

Genesis AI Technology

- Personal health monitoring technology with embedded sensors and artificial intelligence
- Compatible with StarLink accessories

	L	M	P	Matrices: L, M, P Battery Size: 312
Measurement	ANSI/IEC 2cc Coupler	ANSI/IEC 2cc Coupler	ANSI/IEC 2cc Coupler	
Peak OSPL90 (dB SPL)	107	117	120	
HFA OSPL90 (dB SPL)	101	113	116	
Peak Gain (dB)	43	49	60	
HFA Full-On Gain (dB)	36	47	56	
Frequency Range (Hz)	<100-9400	<100-9400	<100-9400	
HFA Frequencies (kHz)	1.0, 1.6, 2.5	N/A	1.0, 1.6, 2.5	
Reference Test Gain (dB)	24	36	39	
Equivalent Input Noise (dB)	26	26	26	
Harmonic Distortion				
500 Hz (%)	<3	<3	<3	
800 Hz (%)	<3	<3	<3	
1600 Hz (%)	<3	<3	<3	
Induction Coil Sensitivity				
HFA SPLITS (ANSI) (dB SPL)	N/A	N/A	N/A	
MASL (IEC) (dB SPL)	N/A	N/A	N/A	
Estimated Battery Life for 16-Hour Day				
Without streaming (days)	Up to 6	Up to 6	Up to 6	
With streaming (days)	Up to 5	Up to 5	Up to 5	
Battery Current (mA)	1.5	1.6	1.6	
Tinnitus Therapy Stimulus				
Max RMS Output (dB SPL)	87	87	87	
Weighted RMS Output Level (dB SPL)	87	87	87	
Max 1/3 Octave Output (dB SPL)	87	87	87	

L

M

P




Starkey | Technical Data

RIC CUSTOM Cased

Receiver-In-Canal

Fitting Range



- L
- M
- P
- UP



Genesis AI 24 | 20 | 16 | 12

Color Guide

Standard Colors



Patient Features

- IP Rating: 68
- Tinnitus Technology
- Wireless Connectivity
- CROS System
- Dual Radio (2.4 GHz + NFMI)

Genesis AI Technology

- Personal health monitoring technology with embedded sensors and artificial intelligence
- Compatible with StarLink accessories

	L	M	P	UP	Matrices: L, M, P, UP Battery Size: 312
Measurement	ANSI/IEC 2cc Coupler	ANSI/IEC 2cc Coupler	ANSI/IEC 2cc Coupler	ANSI/IEC 2cc Coupler	
Peak OSPL90 (dB SPL)	107	117	120	125	
HFA OSPL90 (dB SPL)	101	113	116	122	
Peak Gain (dB)	43	49	60	67	
HFA Full-On Gain (dB)	36	47	56	64	
Frequency Range (Hz)	<100-9400	<100-9400	<100-9400	<100-5800	
HFA Frequencies (kHz)	1.0, 1.6, 2.5	N/A	1.0, 1.6, 2.5	1.0, 1.6, 2.5	
Reference Test Gain (dB)	24	36	39	45	
Equivalent Input Noise (dB)	26	26	26	26	
Harmonic Distortion					
500 Hz (%)	<3	<3	<3	<3	
800 Hz (%)	<3	<3	<3	<3	
1600 Hz (%)	<3	<3	<3	<3	
Induction Coil Sensitivity					
HFA SPLITS (ANSI) (dB SPL)	N/A	N/A	N/A	N/A	
MASL (IEC) (dB SPL)	N/A	N/A	N/A	N/A	
Estimated Battery Life for 16-Hour Day					
Without streaming (days)	Up to 6	Up to 6	Up to 6	Up to 6	
With streaming (days)	Up to 5	Up to 5	Up to 5	Up to 5	
Battery Current (mA)	1.5	1.6	1.6	1.7	
Tinnitus Therapy Stimulus					
Max RMS Output (dB SPL)	87	87	87	87	
Weighted RMS Output Level (dB SPL)	87	87	87	87	
Max 1/3 Octave Output (dB SPL)	87	87	87	87	



Starkey | Technical Data

RIC

Receiver-In-Canal










Starkey | Technical Data

RIC CUSTOM CASED

Receiver-In-Canal

Genesis AI 24 | 20 | 16 | 12



Latency (ms) 4.3

Attack Time (ms) 1

Release Time (ms) 160

Starkey, Genesis, and StarLink are trademarks of Starkey Laboratories, Inc.

©2023 Starkey Laboratories, Inc. All Rights Reserved. 8/23 SPEC3068-01-EN-ST




Starkey | Technical Data

RIC CUSTOM Cased

Receiver-In-Canal

Fitting Range

● UP-



Genesis AI 24 | 20 | 16 | 12

Color Guide

Standard Colors



Patient Features



- IP Rating: 68
- Tinnitus Technology
- Wireless Connectivity
- CROS System
- Dual Radio (2.4 GHz + NFMI)

Genesis AI Technology

- Personal health monitoring technology with embedded sensors and artificial intelligence
- Compatible with StarLink accessories

	UP-	Matrices: UP- Battery Size: 312
--	-----	------------------------------------

Measurement	ANSI/IEC 2cc Coupler
Peak OSPL90 (dB SPL)	119
HFA OSPL90 (dB SPL)	116
Peak Gain (dB)	64
HFA Full-On Gain (dB)	60
Frequency Range (Hz)	<100-5800
HFA Frequencies (kHz)	N/A
Reference Test Gain (dB)	39
Equivalent Input Noise (dB)	26
Harmonic Distortion	
500 Hz (%)	<3
800 Hz (%)	<3
1600 Hz (%)	<3
Induction Coil Sensitivity	
HFA SPLITS (ANSI) (dB SPL)	N/A
MASL (IEC) (dB SPL)	N/A
Estimated Battery Life for 16-Hour Day	
Without streaming (days)	Up to 6
With streaming (days)	Up to 5
Battery Current (mA)	1.6
Tinnitus Therapy Stimulus	
Max RMS Output (dB SPL)	87
Weighted RMS Output Level (dB SPL)	87
Max 1/3 Octave Output (dB SPL)	87



Latency (ms)	4.3
Attack Time (ms)	1
Release Time (ms)	160